

HOW TO PRINT WITH CELLINK A

Overview: Alginate is a naturally derived polysaccharide and has been utilized for diverse applications. This bioink is biocompatible, easy to use, and customizable. CELLINK A can be crosslinked using our crosslinking solution, CaCl₂.

Materials:

CELLINK bioprinter

CELLINK A cartridges

Pluronics or other sacrificial bioink cartridges (also can use CELLINK Support or CELLINK Start)

Conical nozzles

Printing surface, preferably positively charged

Protocol:

1. Mix a high concentration of cells into the CELLINK A solution using CELLINK's Cell Mixer. The Cell Mixer requires about 10 ml of CELLINK A and 1 ml of a high concentration of cells, creating a 10:1 ratio.
2. Attach both syringes to the cell mixer, and attach a 3 cc cartridge to the other end of the Cell Mixer with the piston loaded. Push down on the handle of the Cell Mixer to homogenously mix the two solutions into the 3 cc cartridge.
3. Screw on a sterile nozzle, we recommend using a 25 or 27 G nozzle.
4. Since CELLINK A prints a little runny, print with a sacrificial material in the second print head to help build up the CELLINK A structure. An example code is attached.
5. Once the print has finished, apply enough droplets of 100 µM calcium chloride solution to cover the construct. Crosslinking takes about 5 minutes. Wash with PBS to clear the crosslinking solution.

G-codes:

The following code uses both print heads, and 27G nozzles for both the alginate and the pluronics.

V3_alginateinsidepluronics.gcode

Further Information:

A variation of the CELLINK A bioink is CELLINK A RGD. CELLINK A-RGD Bioink is composed of high purity grade sodium alginate which is bio-functionalized with RGD motifs. CELLINK A-RGD is a

biodegradable bioink and it offers excellent biocompatibility, easy handling, and works universally with a wide range of human cells. CELLINK A-RGD is simply cross-linked with divalent cations such as calcium. You can utilize any of our sacrificial bioinks such as CELLINK Support or CELLINK start when bioprinting with CELLINK A-RGD to create porous structures.

CELLINK A Application Note

References:

N/A

CELLINK AB
Arvid Wallgrens Backe 20
SE 413 46 Gothenburg
Sweden
Phone +46 732 67 00 00

CELLINK LLC
675 W Kendall St.
Cambridge, MA 02142
USA
Phone +1 650 515 5566