

BIOPRINTING PROTOCOL FOR BLOOD VESSELS WITH CELLS

Overview: This protocol is a specific way to create blood vessels using CELLINK bioink and primary human skeletal muscle cells.

Materials:

Autodesk® 123D Design software

Slic3r Software (v1.2.9)

CELLINK bioink

Human smooth muscle cells, SMCs, at a concentration of 10×10^6 cells/ml

3D Cell Culture Media

CELLMIXER Kit

INKREDIBLE 3D Bioprinter by CELLINK

[Conical tip, 250µm ID](#)

CaCl₂ Crosslinking Solution

Protocol:

1. The first step is to design a blue print for the structure. Using Autodesk® 123D Design software, create the 3D model of a tube geometry with the dimensions of 6mm x 10 mm (ODxH)
2. Using Slic3r (v1.2.9), convert the 3D model to a bioprinting protocol and toolpath with the following parameters:
 - Layer height = 0.30mm
 - External perimeters extrusion width = 0.30mm
 - Printing speed, F = 600mm/min

Upload the bioprinting protocol with the following name: "*Blood Vessel_D6xH10_LH_F600.gcode*"

3. The human smooth muscle cells are originally at a concentration of 10×10^6 cell/ml. First, prepare a cell suspension of 33×10^6 SMCs in a volume of 300µL of culture media. See

reference for culture media composition.

<http://dx.doi.org/10.1016/j.biomaterials.2005.10.026>

4. This cell suspension will then be mixed with 3mL of CELLINK^{1,2} bioink using the [CELLMIXER](#) to obtain a final concentration of 10×10^6 cells/ml. This gives a final ratio of 10:1 (bioink to cell suspension). Please watch the video in this link for a detailed illustration on how to do the mixing process:
<https://www.youtube.com/watch?v=CmSYL1-oltI>

5. The following bioprinting parameters can be used with the INKREDIBLE 3D Bioprinter by CELLINK using the pneumatic-driven micro-extrusion technology.
 - [Conical tip, 250µm ID](#) (25 G)
 - Pressure: 15-17 kPa
 - Printing speed: 600 mm/min
 - Printhead temperature: Room temperature (22°C)
 - Printbed temperature: Room temperature (22°C)

6. After the bioprinting process, the cell-laden constructs are crosslinked by submerging in an ionic solution of 100mM CaCl₂ for 6 minutes. The constructs are then rinsed with culture media and incubated in 3D cell culture medium in standard culture conditions (37°C, 5% CO₂ and 95% relative humidity).

7. Bioprinting metrics
 - a. Time for bioprinting: 2 minutes per construct
 - b. Dimensions of bioprinted construct: OD: *6.1mm*, H: *10mm*

8. Post-bioprinting, incubate the cell-laden constructs in 3D cell culture medium and standard culture conditions (37°C, 5% CO₂ and 95% relative humidity) for at least 7 days to analyze the cell viability and morphology.

9. The stability of constructs *in vitro* using macroscopic evaluation showed no change in shape after 7 days of *in vitro* culture.

G-codes:

Blood Vessel_D6xH10_LH03_F600.gcode

Further Information:

blood vessel.stl

CELLINK AB
Arvid Wallgrens Backe 20
SE 413 46 Gothenburg
Sweden
Phone +46 732 67 00 00

CELLINK LLC
675 W Kendall St.
Cambridge, MA 02142
USA
Phone +1 650 515 5566

References:

1. 3D Bioprinting Human Chondrocytes with Nanocellulose–Alginate Bioink for Cartilage Tissue Engineering Applications. Kajsa Markstedt, Athanasios Mantas, Ivan Tournier, Héctor Martínez Ávila, Daniel Hägg, and Paul Gatenholm. *Biomacromolecules* **2015** 16 (5), 1489-1496.
[dx.doi.org/10.1021/acs.biomac.5b00188](https://doi.org/10.1021/acs.biomac.5b00188)
2. 3D bioprinting of human chondrocyte-laden nanocellulose hydrogels for patient-specific auricular cartilage regeneration. Héctor Martínez Ávila, Silke Schwarz, Nicole Rotter, Paul Gatenholm. *Bioprinting* **2016**, Volumes 1–2, 22-35.
[dx.doi.org/10.1016/j.bprint.2016.08.003](https://doi.org/10.1016/j.bprint.2016.08.003)