

Specification Sheet

Crosslinking Agent

Product Description Solution for crosslinking of bioinks, containing 50 mM calcium chloride, **Research Grade**. For research use ONLY. Not intended for *in vitro* diagnostics and *in vivo* uses. Not intended for administration in humans or animals.

Produced under sterile and aseptic conditions.

Batch Number	CL-10100	Expiration date	6 months
Storage	Store at 4-25°C.		
Handling	DO NOT FREEZE. Avoid temperature fluctuations.		
Safety	Handle in accordance with good hygiene and laboratory safety practices.		

Property	Specification	Method
<i>Appearance</i>	Transparent liquid	Visual inspection.
<i>Sterility</i>	Sterile	Tested for the presence of bacteria, fungi and yeast.
<i>Endotoxin level</i>	<10 EU/mL	Limulus Amoebocyte Lysate assay, Pharmacopoeia 2.6.14 "Bacterial endotoxins": Method D, accredited by SWEDAC. Accreditation Certification 1240: ISO 15189, 2010-11-22.
<i>pH</i>	6.5-7.4	Assessed with pH meter.